

EMS Vacuum Immobilisation Price list 2021





RedVac EMS Vacuum Immobilisation

Content	Page	
RedVac made in Austria	3	
The RedVac System	4-5	
Mattresses	6-10	
Head Immobilisation	11	
Extremity Splints	12-16	
Extrication Device	17	
Pediatric	18-19	
Accessories	20	
Pumps & Bags	21	
Thermo Rescue Sack	22	The state of the s
Helicopter Rescue Sack	23	

Developed and manufactured in Austria RedVac



Kohlbrat & Bunz products are developed and manufactured 900 meters above sea level in the middle of the Austrian Alps in its own facilities by its experienced and meticulously working staff.

People living in the Alps know very well that products used under the extreme conditions of alpine rescue have to work perfectly.

To meet these rigorous demands our products are developed in close cooperation with the mountain rescue teams and other demanding emergency services.

This cooperation has brought about a selection of the best materials and optimal product design leading to products of outstanding performance and quality. To insure the best possible quality Kohlbrat & Bunz work in accordance with EN ISO 9001 & EN ISO 13485.



We are very proud of the confidence shown to us and our products by many well-known and demanding emergency services worldwide and are happy to show some examples of their beautiful logos that they have asked us to put on the RedVac products we have made for them.

























The RedVac System

for

PERFECT HYGIENE & EASE OF USE



The unique CLICK-System makes removal of straps and handles for repositioning, replacement, cleaning and disinfection very easy!







The RedVac Multi-Chamber-Systems provides unmatched EASE OF USE

The Multi-Chamber-Systems keeps the internal beads evenly distributed eliminating the need for difficult and time consuming manual spreading of the beads. This enables fast, safe and easy application in all situations.

The RedVac Multi-Chamber-System!

The internal sack system keeps the beads evenly distributed in up to 120 chambers.



The layers of the outer cover are welded together dot by dot creating chambers between the dots keeping the internal beads in their correct position, evenly distributed.







 Optional attachment systems for straps and handles

Cleats







Hook & Loop



Red Stripes



Robust & durable exterior Material

To achieve outstanding durability, abrasion resistance, tightness, easy cleaning and disinfection, comfort and safety for the patient special materials were chosen. It makes it possible to use the RedVac products in a wide temperature range, from -40°C up to +80°C.



 Restraint Straps, made of washable textile or wipe-clean coated fabric

Colour coded patient restraint straps, adjustable from both sides, easy to clean, reposition or replace.



Vacuum Valve, screw-on replaceable
 The valve is robust and easy to use, even when wearing

gloves. It closes automatically when suction stops preventing loss of vacuum if slipping with the pump.

Comes with suitable adapter for most vacuum pumps.



Vacuum Valve, welded-on non-replaceable

The valve is robust and easy to use, even when wearing gloves. It closes automatically when suction stops preventing loss of vacuum if slipping with the pump.



Label System

An indestructible label with individual serial number and date of manufacture is welded onto every product.

Optional labels with property indicator and distributor logo are available.





• X-ray, CT and MRI compatible



CE

All RedVac products are CE marked and conform with the European directive 93/42/EEC for medical products class 1.



ISO9001 - ISO13485 certified by



RedVac EMS Vacuum-Mattresses



Body shape



VM01X02 Cleats with webbing restraint straps **CLICK** with webbing restraint straps VM01CX02 **CLICK** with wipe-clean restraint straps VM01CHX02

Width head/mid/foot: 40-88-50cm Length: 200cm Weight: 7kg

Carry handles: 8 Patient restraint straps: 4

Body shape



VM03X03 Cleats with webbing restraint straps VM03CX03 **CLICK** with webbing restraint straps **CLICK** with wipe-clean restraint straps VM03CHX03

Width head/mid/foot: 50-88-50cm Length: 200cm Weight: 7kg

Carry handles: 8 Patient restraint straps: 4

Body shape. Slits for gurney safety belts in neck area.



Cleats with webbing restraint straps VM03X05 **CLICK** with webbing restraint straps VM03CX05 **CLICK** with wipe-clean restraint straps VM03CHX05

Width head/mid/foot: 50-98-50cm Length: 200cm Weight: 7kg

Carry handles: 8 Patient restraint straps: 4

Body shape. Slits for gurney safety belts in neck area. Glass fiber rod for cervical protection



VM03X07 Cleats with webbing restraint straps VM03CX07 **CLICK** with webbing restraint straps VM03CHX07 **CLICK** with wipe-clean restraint straps

Width head/mid/foot: 50-88-50cm Length: 200cm Weight: 7kg

Carry handles: 8 Patient restraint straps: 4

Body shape - compact size for helicopter tramsportation



VM08X06 Cleats with webbing restraint straps VM08CX06 **CLICK** with webbing restraint straps VM08CHX06 **CLICK** with wipe-clean restraint straps

Width head/mid/foot: 40-78-40cm Length: 195cm Weight: 6kg Carry handles: 8 Patient restraint straps: 4

Cleats with webbing restraint straps VM09X01

CLICK with webbing restraint straps VM09CX01 **CLICK** with wipe-clean restraint straps VM09CHX01 Width head/mid/foot: 40-60-40cm Length: 195cm Weight: 5kg

Carry handles: 8 Patient restraint straps: 4

Body shape with head fixation.

Slits for gurney safety belts in neck area. Glass fiber rod for cervical protection.



Cleats with webbing restraint straps VM200X06 **CLICK** with webbing restraint straps VM200CX06 **CLICK** with wipe-clean restraint straps VM200CHX06

Width head/mid/foot: 40-88-60cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 6

VM201X10 Cleats with webbing restraint straps **CLICK** with webbing restraint straps VM201CX10 VM201CHX10 **CLICK** with wipe-clean restraint straps

Width head/mid/foot: 40-100-60cm Length: 210cm Weight: 8kg Carry handles: 8 Patient restraint straps: 6



Head fixation VM500X03 Cleats with webbing restraint straps VM500CX03 **CLICK** with webbing restraint straps CLICK with wipe-clean restraint straps VM500CHX03 Width head/mid/foot: 60-88-88cm Length: 200cm Weight: 7,5kg Carry handles: 8 Patient restraint straps: 5 VM501X03 Cleats with webbing restraint straps VM501CX03 **CLICK** with webbing restraint straps **CLICK** with wipe-clean restraint straps VM501CHX03 Width head/mid/foot: 60-100-100cm Length: 210cm Weight: 8,5kg Carry handles: 8 Patient restraint straps: 5 Cleats with webbing restraint straps VM502X01 **CLICK** with webbing restraint straps VM502CX01 VM502CHX01. **CLICK** with wipe-clean restraint straps

Head fixation and slits in neck area for gurney safety belts. Glassfiber rod for cervical protection.

Carry handles: 8 Patient restraint straps: 5

Width head/mid/foot: 60-80-80cm Length: 200cm Weight: 7kg

ricaa rixatiori aria siits	III HEEK area for garri	cy surety bents. Glassiber roa for cervical protection.
	VM5000X01	Cleats with webbing restraint straps
	VM5000CX01 VM5000CHX01	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps
		65-88-88cm Length: 200cm Weight: 7kg ient restraint straps: 5
A S P P P	VM5001X01	Cleats with webbing restraint straps
September 18	VM5001CX01 VM5001CHX01	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps
		65-100-100cm Length: 210cm Weight: 8kg attent restraint straps: 5
	VM5002X01	Cleats with webbing restraint straps
	VM5002CX01 VM5002CHX01	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps
		65-80-80cm Length: 200cm Weight: 6kg ient restraint straps: 5

Body shape with head fixation. Slits for gurney safety belts. Glass fiber rod for cervical protection. Additional outer protective layer with 10 integrated handles.

riadicional dator proto-	cuve layer with 10 integrated namines.
	VM7000X03 Cleats with webbing restraint straps
n o	VM7000CX03 CLICK with webbing restraint straps VM7000CHX03 CLICK with wipe-clean restraint straps
CIL	Width: 130cm Length: 205cm Weight: 10kg Carry handles non-removable: 10 Carry handles removable: 4 Patient restraint straps: 5
	VM7000X04 Cleats with webbing restraint straps
	VM7000CX04 CLICK with webbing restraint straps VM7000CHX04 CLICK with wipe-clean restraint straps
	Width: 130cm Length: 205cm Weight: 10kg Carry handles non-removable: 10 Patient restraint straps: 5
\n + 1 1	VM7001X01 Cleats with webbing restraint straps
-	VM7001CX01 CLICK with webbing restraint straps VM7001CHX01 CLICK with wipe-clean restraint straps
HIT	Width: 108cm Length: 205cm Weight: 9,5kg Carry handles non-removable: 10 Carry handles removable: 4 Patient restraint straps: 5
	VM7001X04 Cleats with webbing restraint straps
	VM7001CX04 CLICK with webbing restraint straps VM7001CHX04 CLICK with wipe-clean restraint straps
THIT	Width: 110cm Length: 205cm Weight: 9,5kg Carry handles non-removable: 10 Patient restraint straps: 5



Head fixation. Slits for gurney safety belts in neck area. Glassfiber rod for cervical protection



VM11X04 Cleats with webbing restraint straps
VM11CX04 CLICK with webbing restraint straps
VM11CHX04 CLICK with wipe-clean restraint straps


Width head/mid/foot: 80-88-80cm Length: 200cm Weight: 7kg

Carry handles: 8 Restraint straps: 4

VM1101X01 Cleats with webbing restraint straps
VM1101CX01 CLICK with webbing restraint straps
VM1101CHX04 CLICK with wipe-clean restraint straps

Width head/mid/foot: 80-100-80cm Length: 210cm Weight: 8kg

Carry handles: 12 Restraint straps: 4



VM1111X01 Cleats with webbing restraint straps
VM1111CX01 CLICK with webbing restraint straps
VM1111CHX01 CLICK with wipe-clean restraint straps

Width head/mid/foot: 80-108-80cm Length: 210cm Weight: 8,5kg

Carry handles: 12 Restraint straps: 4

VM1112X01 Cleats with webbing restraint straps
VM1112CX01 CLICK with webbing restraint straps
VM1112CHX01 CLICK with wipe-clean restraint straps

Width head/mid/foot: 80-108-80cm Length: 220cm Weight: 8,5kg

Carry handles: 8 Restraint straps: 4

Rectangular shape. Slits for gurney safety belts in the neck area Glassfiber rod for cervical protection.

VM11620X01	Cleats with webbing restraint straps
VM11620CX01 VM11620CHX01	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps
VWITTOZOCITACI	CLICK with wipe-clean restraint straps

Width: 60cm Length: 200cm Weight: 5kg

Carry handles: 8 Restraint straps: 4

VM11720X01 Cleats with webbing restraint straps
VM11720CX01 CLICK with webbing restraint straps
VM11720CHX01 CLICK with wipe-clean restraint straps

Width: 70cm Length: 200cm Weight: 6kg

Carry handles: 8 Restraint straps: 4



VM11820X01 Cleats with webbing restraint straps
VM11820CX01 CLICK with webbing restraint straps
VM11820CHX01 CLICK with wipe-clean restraint straps

Width: 80cm Length: 200cm Weight: 7kg

Carry handles: 8 Restraint straps: 4

VM11920X01 Cleats with webbing restraint straps
VM11920CX01 CLICK with webbing restraint straps
VM11920CHX01 CLICK with wipe-clean restraint straps

Width: 90cm Length: 200cm Weight: 7,5kg

Carry handles: 8 Restraint straps: 4

VM1201X01 Cleats with webbing restraint straps
VM1201CX01 CLICK with webbing restraint straps
VM1201CHX01 CLICK with wipe-clean restraint straps

Width: 100cm Length: 210cm Weight: 8,5kg

Carry handles: 8 Restraint straps: 4



Tapering shape at hea	nd and foot area		
12.12.1.2.18	VM10X02	Cleats with webbing restraint straps	
	VM10CX02 VM10CHX02	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps	
3 4 2 16	Width head/mid/foot: 80 Carry handles: 8 Restra	0-88-80cm Length: 200cm Weight: 7kg aint straps: 4	
	VM8021X01	Cleats with webbing restraint straps	
8. 7. 9. 9	VM8021CX01 VM8021CHX01	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps	
Rectangular shape	Width head/mid/foot: 80 Carry handles: 12 Rest	0-100-80cm Length: 210cm Weight: 8,5kg raint straps: 4	
	VM8121X01	Cleats with webbing restraint straps	
	VM8121CX01 VM8121CHX01	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps	
	Width head/mid/foot: 80 Carry handles: 12 Rest	0-108-80cm Length: 210cm Weight: 9kg raint straps: 4	
	l		
Rectangular shape	,		
	VM620X01	Cleats with webbing restraint straps	
	VM620CX01 VM620CHX01	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps	
	Width: 60cm Lengtl Carry handles: 8 Restra	n: 200cm Weight: 6kg aint straps: 4	
	VM720X01	Cleats with webbing restraint straps	
20 8 8 8 8	VM720CX01 VM720CHX01	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps	
	Width: 70cm Lengt Carry handles: 8 Restra	h: 200cm Weight: 6,5kg aint straps: 4	
13141414	VM820X01	Cleats with webbing restraint straps	
	VM820CX01 VM820CHX01	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps	
	Width: 80cm Lengt Carry handles: 8 Restra	h: 200cm Weight: 7kg aint straps: 4	
	VM1021X01	Cleats with webbing restraint straps	
	VM1021CX01 VM1021CHX01	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps	
	Width: 100cm Lengt Carry handles: 12 Rest	h: 210cm Weight: 8,5kg raint straps: 4	
	VM1121X01	Cleats with webbing restraint straps	
	VM1121CX01 VM1121CHX01	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps	
	Milelthe 400 am	h. 240am Wainht. Oka Handlas 42	l

Weight: 9kg

Handles: 12


Width: 108cm Length: 210cm Carry handles: 12 Restraint straps: 4



For obese patients Slits for gurney safety be	elts in the neck are	ea Glassfiber rod for cervical protection.	
	VM1421X01	Cleats with webbing restraint straps	
4 6 6 6	VM1421CX01 VM1421CHX01	CLICK with webbing restraint straps CLICK with wipe-clean restraint straps	

The RedVac Combi-Vacuum-Mattress

2 in 1 - Combines the standard mattress and the vacuum mattress





The RedVac Vacuum-Head-Immobiliser

- for the Spine Board & the Scoop Stretcher

- Quick and easy positioning of the immobiliser
- No manipulation of head and cervical spine
- Adjusts automatically to any head shape and size
- Ears remain free for unrestricted hearing and medical examination
- Solid and comfortable immobilisation



Specification		
	Head Immobiliser VK98X07 Width: 38cm Length: 31cm Weight: 0,7kg Base plates GP08 for scoop stretcher	
	GP11 for back board	
	Head Immobiliser with buckles for attachment VK94X17	
	Width: 38cm Length: 31cm Weight: 0,7kg	



RedVac Vacuum Extrication Device

- the answer to your extrication challenges!
- fits every patient, pregnant – overweight … no problem!
- comfortable no respiratory impairment!

Specification	Price in
VED24X02	
Width, head: 40cm Width, middle: 106cm Length: 94cm Weight: 2,5kg	



RedVac EMS Vacuum Extremity Splints



Optional restraint systems:

CLICK-System Hook & Loop Red Stripes Cleats









Patient straps:
Webbing of PP
Wipe-clean of PVC
Choice of colours



Angulate design		
	Forearm Splints VS60X01 Red Stripes with Webbing 1-Chamber VS61X01 Red Stripes with Webbing Multi-Chamber-Sack Restraint straps: 3 Width: 35cm Length: 50cm Weight: 0,33kg	
	Arm Splints	
	VS62X01 Red Stripes with Webbing 1-Chamber VS63X03 Red Stripes with Webbing Multi-Chamber-Sack	
	Restraint straps: 5 Width top/foot: 43-70cm Length: 56cm Weight: 0,9kg	
1	Leg Splints VS64X21 Red Stripes with Webbing 1-Chamber VS65X21 Red Stripes with Webbing Multi-Chamber-Sack Restraint straps: 5 Width top/foot: 64-70cm Length: 103cm Weight: 1,9kg	
The state of		
NIM 24	Forearm Splints VS60CX05 CLICK with Webbing VS60CHX05 CLICK with Wipe-Clean Straps VS61CX05 CLICK with Webbing Multi-Chamber-Sack	
44-7	VS61CHX05 CLICK with Wipe-Clean Straps Multi-Chamber-Sack VS61CPX05 CLICK with Webbing DOT-Chamber VS61CPHX05 CLICK with Wipe-Clean Straps DOT-Chamber	
	Restraint straps: 3 Width: 35cm Length: 50cm Weight: 0,33kg	
1:1	Arm Splints VS62CX05 CLICK with Webbing 1-Chamber VS62CHX05 CLICK with Wipe-Clean Straps 1-Chamber	
	VS63CX05 CLICK with Webbing Multi-Chamber-Sack VS63CHX05 CLICK with Wipe-Clean Straps Multi-Chamber-Sack VS63CPX05 CLICK with Webbing DOT-Chamber	
	VS63CPHX05 CLICK with Wipe-Clean Straps DOT-Chamber	
-	Restraint straps: 5 Width top/foot: 43-70cm Length: 56cm Weight: 0,9kg	
	Leg Splints VS64CX05 CLICK with Webbing 1-Chamber VS64CHX05 CLICK with Wipe-Clean Straps 1-Chamber	
	VS65CX23 CLICK with Webbing Multi-Chamber-Sack VS65CHX23 CLICK with Wipe-Clean Straps Multi-Chamber-Sack	
	VS65CPX23 CLICK with Webbing DOT-Chamber VS65CPHX23 CLICK with Wipe-Clean Straps DOT-Chamber	
	Restraint straps: 5 Width top/foot: 64-70cm Length: 103cm Weight: 1,9kg	
00	Forearm Splints VS60X20V Hook & Loop straps 1-Chamber VS61X20V Hook & Loop straps Multi-Chamber-Sack VS61PX20V Hook & Loop straps DOT-Chamber	
-	Restraint straps: 2 Width: 35cm Length: 50cm Weight: 0,33kg	
	Arm Splints VS62X20V Hook & Loop straps 1-Chamber VS63X20V Hook & Loop straps Multi-Chamber-Sack VS63PX20V Hook & Loop straps DOT-Chamber	
	Restraint straps: 4 Width top/foot: 43-70cm Length: 56cm Weight: 0,9kg Welded-on non-removable valve	
\t	Leg Splints VS64X20V Hook & Loop straps 1-Chamber VS65X20V Hook & Loop straps Multi-Chamber-Sack VS65PX20V Hook & Loop straps DOT-Chamber	
	Restraint straps: 5 Width top/foot: 64-70cm Length: 103cm Weight: 1,9kg Welded-on non-removable valve	



	Arm splints	
N market	VA14X02 Red Stripes with Webbing 1-Chamber VA14X01 Red Stripes with Webbing Multi-Chamber-Sack	
	Restraint straps: 4 Width: top/hand: 50-30cm Length: 67cm Weight: 0,9kg	
	Ankle splints	
	VK15X02 Red Stripes with Webbing 1-Chamber VK15X01 Red Stripes with Webbing Multi-Chamber-Sack	
Control Control	Restraint straps: 4 Width: 47cm Length: 57cm Weight: 0,9kg	
	Leg splints	
N. Comments	VL21X02 Red Stripes with Webbing 1-Chamber VL21X01 Red Stripes with Webbing Multi-Chamber-Sack	
	Restraint straps: 5 Width: top/hand: 70-40cm Length: 110cm Weight: 1,9kg	
	Large leg splints	
	VL18X02 Red Stripes with Webbing 1-Chamber VL18X01 Red Stripes with Webbing Multi-Chamber-Sack	
And the second	VL18X01 Red Stripes with Webbing Multi-Chamber-Sack Restraint straps: 5 Width: top/hand: 70-45cm Length: 135cm Weight: 2,1kg	
	Arm splints	
	VA14X04 CLEAT with webbing 1-Chamber	
	VA14X03 CLEAT with Webbing Multi-Chamber-Sack	
	VA14CX04 CLICK with webbing 1-Chamber VA14CHX04 CLICK with Wipe-Clean Straps 1-Chamber	
BER 1881	VA14CX03 CLICK with Webbing Multi-Chamber-Sack	
4 7 7	VA14CHX03 CLICK with Wipe-Clean Straps Multi-Chamber-Sack VA14CPX03 CLICK with Webbing DOT-Chamber	
	VA14CPX03 CLICK with Webbing DOT-Chamber VA14CPHX03 CLICK with Wipe-Clean Straps DOT-Chamber	
	Restraint straps: 4 Width: top/hand: 50-30cm Length: 67cm Weight: 0,9kg	
	Ankle splints	
	VK15X04 CLEAT with webbing 1-Chamber VK15X03 CLEAT with Webbing Multi-Chamber-Sack	
	VK15X03 CLEAT with Webbing Multi-Chamber-Sack VK15CX04 CLICK with webbing 1-Chamber	
*	VK15CHX04 CLICK with Wipe-Clean Straps 1-Chamber	
	VK15CX03 CLICK with Webbing Multi-Chamber-Sack VK15CHX03 CLICK with Wipe-Clean Straps Multi-Chamber-Sack	
	VK15CPX03 CLICK with Webbing DOT-Chamber	
	VK15CPHX03 CLICK with Wipe-Clean Straps DOT-Chamber	
	Restraint straps: 4 Width: 47cm Length: 57cm Weight: 0,9kg	
11	Leg splints VL21X04 CLEAT with webbing 1-Chamber	
Mar. 100	VL21X03 CLEAT with Webbing I-Chamber VL21X03 CLEAT with Webbing Multi-Chamber-Sack	
	VL21CX04 CLICK with webbing 1-Chamber VL21CHX04 CLICK with Wipe-Clean Straps 1-Chamber	
	VL21CX03 CLICK with Webbing 1-Chamber Sack	
	VL21CHX03 CLICK with Wipe-Clean Straps Multi-Chamber-Sack	
	VL21CPX03 CLICK with Webbing DOT-Chamber VL21CPHX03 CLICK with Wipe-Clean Straps DOT-Chamber	
	Restraint straps: 5 Width: top/hand: 70-40cm Length: 110cm Weight: 1,9kg	
	Large leg splints VL18X04 CLEAT with webbing 1-Chamber	
	VL18X03 CLEAT with Webbing Multi-Chamber-Sack	
11	VL18CX04 CLICK with webbing 1-Chamber VL18CHX04 CLICK with Wipe-Clean Straps 1-Chamber	
A STATE OF THE PERSON NAMED IN	VL18CX03 CLICK with Webbing Multi-Chamber-Sack	
Transaction of the last of the	VL18CHX03 CLICK with Wipe-Clean Straps Multi-Chamber-Sack	
	VL18CPX03 CLICK with Webbing DOT-Chamber VL18CPHX03 CLICK with Wipe-Clean Straps DOT-Chamber	
	Restraint straps: 5 Width: top/hand: 70-45cm Length: 135cm Weight: 2,1kg	
	25.18g 1.000m 10.00m 25.18g	



	Forearm Splints	
	VS100CKX02 Loop strap with CLICK 1-Chamber VS100CKX01 Loop strap with CLICK Multi-Chamber-Sack VS100CPKX01 Loop strap with CLICK DOT-Chamber	
	Restraint straps: 2 Width: 35cm Length: 53cm Weight: 0,5kg	
	Arm splints	
	VS120CKX02 Loop strap with CLICK 1-Chamber VS120CKX01 Loop strap with CLICK Multi-Chamber-Sack VS120CPKX01 Loop strap with CLICK DOT-Chamber Restraint straps: 4 Width: 63-47cm Length: 74cm Weight: 1,1kg	
	Leg Splints VS140CKX02 Loop strap with CLICK 1-Chamber VS140CKX01 Loop strap with CLICK Multi-Chamber-Sack VS140CPKX01 Loop strap with CLICK DOT-Chamber	
E	Restraint straps: 6 Width: 70-36cm Length: 110cm Weight: 2,1kg	

	Forearm Splin	ts		
-	VS100CX02 VS100CHX02	Webbing with CLICK Wipe-Clean-Strap with CLICK	1-Chamber 1-Chamber	
	VS100CX01 VS100CHX01	Webbing with CLICK Wipe-Clean-Strap with CLICK	Multi-Chamber-Sack Multi-Chamber-Sack	
-	VS100CPX01 VS100CPHX01	Webbing with CLICK Wipe-Clean-Strap with CLICK	DOT-Chamber DOT-Chamber	
	Restraint straps	: 3 Width: 35cm Length: 53cm We	eight: 0,5kg	
Arm	Arm splints			
- 40	VS120CX02 VS120CHX02	Webbing with CLICK Wipe-Clean-Strap with CLICK	1-Chamber 1-Chamber	
-	VS120CX01 VS120CHX01	Webbing with CLICK Wipe-Clean-Strap with CLICK	Multi-Chamber-Sack Multi-Chamber-Sack	
	VS120CPX01 VS120CPHX01	Webbing with CLICK Wipe-Clean-Strap with CLICK	DOT-Chamber DOT-Chamber	
	Restraint straps	: 4 Width: 63-47cm Length: 74cm	Weight: 1,1kg	
Bein	Leg Splints			
	VS140CX02 VS140CHX02	Webbing with CLICK Wipe-Clean-Strap with CLICK	1-Chamber 1-Chamber	
-	VS140CX01 VS140CHX01	Webbing with CLICK Wipe-Clean-Strap with CLICK	Multi-Chamber-Sack Multi-Chamber-Sack	
1	VS140CPX01 VS140CPHX01	Webbing with CLICK Wipe-Clean-Strap with CLICK	DOT-Chamber DOT-Chamber	
-	Restraint straps	: 6 Width: 70-36cm Length: 110cm	Weight: 2,1kg	



Forearm Splints		
VS100CKX04	Loop strap with CLICK	1-Chamber
VS100CKX03	Loop strap with CLICK	Multi-Chamber-Sack
VS100CPKX03	Loop strap with CLICK	DOT-Chamber
Restraint straps: 2 Wid	dth: 35cm Length: 53cm We	eight: 0,5kg
Arm Splints		
7 VS120CKX04	Loop strap with CLICK	1-Chamber
	Loop strap with CLICK	Multi-Chamber-Sack
VS120CPKX03	Loop strap with CLICK	DOT-Chamber
Restraint straps: 4 Wid	dth: 63-47cm Length: 74cm	Weight: 1,1kg
Leg Splints		
VS140CKX04	Loop strap with CLICK	1-Chamber
MARKS -	• •	Multi-Chamber-Sack
VS140CPKX03	Loop strap with CLICK	DOT-Chamber
The state of the s		
Restraint straps: 6 Wid	ith: 70-36cm Length: 110cm	Weight: 2,1kg
Restraint straps: 6 Wid	dth: 70-36cm Length: 110cm	Weight: 2,1kg
Forearm Splints VS100CX04 Web	bbing with CLICK be-Clean-Strap with CLICK	1-Chamber 1-Chamber
Forearm Splints VS100CX04 Web VS100CHX04 Wip VS100CX03 Web	bbing with CLICK	1-Chamber
Forearm Splints VS100CX04 Web VS100CHX04 Wip VS100CX03 Web VS100CHX03 Wip VS100CPX03 Web	bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK	1-Chamber 1-Chamber Multi-Chamber-Sack
Forearm Splints VS100CX04 Web VS100CHX04 Wip VS100CHX03 Web VS100CHX03 Wip VS100CPX03 Web VS100CPX03 Web VS100CPX03 Web	bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK	1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber
Forearm Splints VS100CX04 Web VS100CHX04 Wip VS100CHX03 Web VS100CHX03 Wip VS100CPX03 Web VS100CPX03 Web VS100CPX03 Web	bbing with CLICK pe-Clean-Strap with CLICK bbing with CLICK pe-Clean-Strap with CLICK bbing with CLICK	1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber
Forearm Splints V\$100CX04 Web V\$100CHX04 Wip V\$100CHX03 Web V\$100CHX03 Wip V\$100CPX03 Web V\$100CPX03 Web V\$100CPHX03 Wip Restraint straps: 3 Wid Arm Splints V\$120CX04 Web	bbing with CLICK pe-Clean-Strap with CLICK bbing with CLICK pe-Clean-Strap with CLICK bbing with CLICK	1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber
Forearm Splints VS100CX04 Web VS100CHX04 Wip VS100CX03 Web VS100CHX03 Wip VS100CPX03 Web VS100CPX03 Wip Restraint straps: 3 Wid Arm Splints VS120CX04 Web VS120CX04 Wip VS120CX03 Web	bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK dth: 35cm Length: 53cm We bbing with CLICK bbing with CLICK	1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber eight: 0,5kg
Forearm Splints V\$100CX04 Web V\$100CHX04 Wip V\$100CHX03 Web V\$100CPX03 Web V\$100CPX03 Web V\$100CPHX03 Wip Restraint straps: 3 Wid Arm Splints V\$120CX04 Web V\$120CHX04 Wip V\$120CHX03 Wip V\$120CHX03 Wip V\$120CHX03 Wip V\$120CHX03 Web V\$120CPX03 Web	bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK dth: 35cm Length: 53cm We bbing with CLICK be-Clean-Strap with CLICK be-Clean-Strap with CLICK bbing with CLICK bbing with CLICK bbing with CLICK	1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber eight: 0,5kg 1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber
Forearm Splints V\$100CX04 Web V\$100CHX04 Wip V\$100CHX03 Web V\$100CPX03 Web V\$100CPX03 Wip Restraint straps: 3 Wid Arm Splints V\$120CX04 Web V\$120CHX04 Wipo	bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK dth: 35cm Length: 53cm We bbing with CLICK be-Clean-Strap with CLICK	1-Chamber 1-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber eight: 0,5kg 1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber
Forearm Splints V\$100CX04 Web V\$100CHX04 Wip V\$100CHX03 Web V\$100CPX03 Web V\$100CPX03 Wip Restraint straps: 3 Wid Arm Splints V\$120CX04 Web V\$120CHX04 Wipo V\$120CHX06 Wipo V\$120CHX06 Wipo	bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK dth: 35cm Length: 53cm We bbing with CLICK be-Clean-Strap with CLICK be-Clean-Strap with CLICK bbing with CLICK bbing with CLICK bbing with CLICK	1-Chamber 1-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber eight: 0,5kg 1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber
Forearm Splints V\$100CX04 Web V\$100CHX04 Wip V\$100CHX03 Web V\$100CPX03 Web V\$100CPX03 Wip Restraint straps: 3 Wid Arm Splints V\$120CX04 Web V\$120CHX04 Wipo	bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK dth: 35cm Length: 53cm We bbing with CLICK be-Clean-Strap with CLICK	1-Chamber 1-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber eight: 0,5kg 1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber
Forearm Splints V\$100CX04 Web V\$100CHX04 Wip V\$100CHX03 Web V\$100CPX03 Web V\$100CPX03 Wip Restraint straps: 3 Wid Arm Splints V\$120CX04 Web V\$120CX04 Web V\$120CHX04 Wip V\$120CHX04 Wip V\$120CHX04 Wip V\$120CHX03 Wip Restraint straps: 4 Wid Leg Splints V\$140CX04 Web	bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK dth: 35cm Length: 53cm We bbing with CLICK be-Clean-Strap with CLICK	1-Chamber 1-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber eight: 0,5kg 1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber
Forearm Splints V\$100CX04 Web V\$100CHX04 Wip V\$100CX03 Web V\$100CPX03 Wip V\$100CPX03 Wip Restraint straps: 3 Wid Arm Splints V\$120CX04 Web V\$120CX04 Wip V\$120CHX04 Wip V\$120CHX04 Wip V\$120CHX03 Wip Restraint straps: 4 Wid Leg Splints V\$140CX04 Web V\$140CX04 Web V\$140CX04 Wip V\$140CX03 Web	bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK dth: 35cm Length: 53cm We bbing with CLICK be-Clean-Strap with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK	1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber eight: 0,5kg 1-Chamber 1-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber UT-Chamber Multi-Chamber Multi-Chamber Multi-Chamber
Forearm Splints V\$100CX04 Web V\$100CHX04 Wip V\$100CHX03 Wip V\$100CPX03 Web V\$100CPX03 Wip Restraint straps: 3 Wid Arm Splints V\$120CX04 Web V\$120CHX04 Wip V\$120CHX04 Wip V\$120CHX03 Wip Restraint straps: 4 Wid Leg Splints V\$120CPX03 Web V\$120CPX03 Wip Restraint straps: 4 Wid Leg Splints V\$140CX04 Web V\$140CHX04 Wip V\$140CHX04 Wip V\$140CHX04 Wip V\$140CHX04 Wip V\$140CHX04 Wip V\$140CHX04 Wip	bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK dth: 35cm Length: 53cm We bbing with CLICK be-Clean-Strap with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK	1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber eight: 0,5kg 1-Chamber 1-Chamber-Sack Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber UT-Chamber DOT-Chamber DOT-Chamber DOT-Chamber
Forearm Splints V\$100CX04 Web V\$100CHX04 Wip V\$100CHX03 Wip V\$100CPX03 Web V\$100CPX03 Wip Restraint straps: 3 Wid Arm Splints V\$120CX04 Web V\$120CHX04 Wip V\$120CHX03 Wip V\$120CHX03 Wip V\$120CPX03 Web V\$120CPX03 Web V\$120CPX03 Wip Restraint straps: 4 Wid Leg Splints V\$140CX04 Web V\$140CX04 Wip V\$140CX04 Wip V\$140CX03 Web	bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bth: 35cm Length: 53cm We bbing with CLICK be-Clean-Strap with CLICK be-Clean-Strap with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK bbing with CLICK be-Clean-Strap with CLICK bbing with CLICK	1-Chamber 1-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber eight: 0,5kg 1-Chamber 1-Chamber-Sack Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber DOT-Chamber Weight: 1,1kg 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack DOT-Chamber OT-Chamber



Splints with Hook &	Loop	
	Arm Splint VS101X01 Hook & Loop straps with Multi-Chamber-Sack Restraint straps: 2 Width: Top/Hand: 40-30cm Length: 73cm Weight: 0,7kg	
	VS103X01 Hook & Loop straps with Multi-Chamber-Sack Restraint straps: 4 Width: Top/Foot: 60-18cm Length: 96cm Weight: 1,9kg	
	VS105X01 Hook & Loop straps with Multi-Chamber-Sack Restraint straps: 5 Width: Top/Foot: 73-18cm Length: 133cm Weight: 2,2kg	

VS101X02 Hook & Loop straps with Multi-Chamber-Sack Restraint straps: 2 Width: Top/Hand: 40-30cm Length: 73cm Weight: 0,7kg	
VS103X02 Hook & Loop straps with Multi-Chamber-Sack Restraint straps: 4 Width: Top/Foot: 60-18cm Length: 96cm Weight: 1,9kg	
VS105X02 Hook & Loop straps with Multi-Chamber-Sack Restraint straps: 5 Width: Top/Foot: 73-18cm Length: 133cm Weight: 2,2kg	





Forearm splint

VS66X03 **Cleats with Webbing** 1-Chamber VS66CX03 **CLICK** with Webbing 1-Chamber VS66CHX03 1-Chamber **CLICK with Wipe-Clean Straps**

VS67X03 Cleats with Webbing Multi-Chamber-Sack VS67CX03 **CLICK with Webbing** Multi-Chamber-Sack VS67CHX03 **CLICK with Wipe-Clean Straps** Multi-Chamber-Sack VS67CPX03 **CLICK** with Webbing **DOT-Chamber** VS67CPHX03 **DOT-Chamber CLICK with Wipe-Clean Straps**

Restraint straps: 3 Width: 30cm Length: 35cm Weight: 0.3kg



Arm splint

VS68X04 Cleats with Webbing 1-Chamber **CLICK** with Webbing VS68CX04 1-Chamber VS68CHX04 **CLICK with Wipe-Clean Straps** 1-Chamber

Cleats with Webbing Multi-Chamber-Sack VS69X03 **CLICK with Webbing** Multi-Chamber-Sack VS69CX03 VS69CHX03 **CLICK with Wipe-Clean Straps** Multi-Chamber-Sack **CLICK** with Webbing **DOT-Chamber** VS69CPX03

CLICK with Wipe-Clean Straps DOT-Chamber VS69CPHX03 Restraint straps: 4 Width: top/hand: 65-44cm Length: 70cm Weight: 1,5kg



Leg splint

VS70X03 Cleats with Webbing 1-Chamber VS70CX03 **CLICK** with Webbing 1-Chamber **VS70CHX03** CLICK with Wipe-Clean Straps 1-Chamber

VS71X03 Cleats with Webbing Multi-Chamber-Sack VS71CX03 CLICK with Webbing Multi-Chamber-Sack **VS71CHX03 CLICK with Wipe-Clean Straps** Multi-Chamber-Sack **VS71CPX03** DOT-Chamber **CLICK** with Webbing VS71CPHX03 **CLICK with Wipe-Clean Straps DOT-Chamber**

Restraint straps: 5 Width: top/foot: 70-55cm Length: 90cm Weight: 1,9kg



Forearm splint

VS66X04 Cleats with Webbing 1-Chamber VS66CX04 1-Chamber **CLICK with Webbing** VS66CHX04 **CLICK with Wipe-Clean Straps** 1-Chamber

VS67X04 **Cleats with Webbing** Multi-Chamber-Sack VS67CX04 **CLICK with Webbing** Multi-Chamber-Sack VS67CHX04 CLICK with Wipe-Clean Straps Multi-Chamber-Sack VS67CPX04 **CLICK** with Webbing **DOT-Chamber** VS67CPHX04 **CLICK with Wipe-Clean Straps DOT-Chamber**

Restraint straps: 3 Width: 30cm Length: 35cm Weight: 0,3kg



Arm splint

VS68X05 Cleats with Webbing 1-Chamber VS68CX05 **CLICK with Webbing** 1-Chamber VS68CHX05 **CLICK with Wipe-Clean Straps** 1-Chamber

VS69X04 Cleats or Stripes with Webbing Multi-Chamber-Sack VS69CX04 **CLICK** with Webbing Multi-Chamber-Sack VS69CHX04 **CLICK with Wipe-Clean Straps** Multi-Chamber-Sack **CLICK** with Webbing VS69CPX04 **DOT-Chamber** VS69CPHX04

DOT-Chamber

Restraint straps: 4 Width: top/hand: 65-44cm Length: 70cm Weight: 1,5kg



Leg splint

VS70X04 Cleats with Webbing 1-Chamber VS70CX04 CLICK with Webbing 1-Chamber VS70CHX04 **CLICK with Wipe-Clean Straps** 1-Chamber

CLICK with Wipe-Clean Straps

VS71X04 Cleats with Webbing Multi-Chamber-Sack VS71CX04 **CLICK with Webbing** Multi-Chamber-Sack VS71CHX04 **CLICK with Wipe-Clean Straps** Multi-Chamber-Sack VS71CPX04 CLICK with Webbing **DOT-Chamber** VS71CPHX04 CLICK with Wipe-Clean Straps **DOT-Chamber**

Restraint straps: 5 Width: top/foot: 70-55cm Length: 90cm Weight: 1,9kg

RedVac Pediatric




Vacuum-Mattresses				
	VM710X11 VM710PX11 4 built in handles	No patient restraint straps 1-Cham No patient restraint straps DOT-Ch Width: 70cm Length: 100cm Wei		
18.18	VM710X09 VM710X08 VM710CX09 VM710CX08 VM710CPX08 VM710CHX09 VM710CHX08 VM710CPHX08 Carry handles: 4 Width: 70cm Le	Cleats with Webbing Cleats with Webbing CLICK with Webbing CLICK with Webbing CLICK with Webbing CLICK with Wipe-Clean Straps CLICK with Wipe-Clean Straps CLICK with Wipe-Clean Straps CLICK with Wipe-Clean Straps Restraint straps: 3 ngth: 100cm Weight: 3,5kg	1-Chamber Multi-Chamber-Sack 1-Chamber Multi-Chamber-Sack DOT-Chamber 1-Chamber Multi-Chamber-Sack DOT-Chamber-Sack	
18. [8]	VM710X14 VM710X13 VM710CX14 VM710CX13 VM710CPX13 VM710CHX14 VM710CHX13 VM710CPHX13 Carry handles: 4 Width: 70cm Le	Cleats with Webbing Cleats with Webbing CLICK with Webbing CLICK with Webbing CLICK with Webbing CLICK with Wipe-Clean Straps CLICK with Wipe-Clean Straps CLICK with Wipe-Clean Straps CLICK with Wipe-Clean Straps Restraint straps: 3 Ingth: 100cm Weight: 4kg	1-Chamber Multi-Chamber-Sack 1-Chamber Multi-Chamber-Sack DOT-Chamber 1-Chamber Multi-Chamber-Sack DOT-Chamber	
18-18	VM712X10 VM712X09 VM712CX10 VM712CX09 VM712CPX09 VM712CHX10 VM712CHX09 VM712CPHX09 Carry handles: 4 Width: 70cm Le	Cleats with Webbing Cleats with Webbing CLICK with Webbing CLICK with Webbing CLICK with Webbing CLICK with Wipe-Clean Straps CLICK with Wipe-Clean Straps CLICK with Wipe-Clean Straps CLICK with Wipe-Clean Straps Restraint straps: 3 ngth: 120cm Weight: 4,5kg	1-Chamber Multi-Chamber-Sack 1-Chamber Multi-Chamber-Sack DOT-Chamber 1-Chamber Multi-Chamber-Sack DOT-Chamber	
Pediatric Vacuum Sp	Arm splint, angled VS56X20V Restraint straps: 3	Hook & Loop straps 0-28cm Length: 42cm Weight: 0,6 kg	1-Chamber	
	Leg Splint VS58X20V Restraint straps: 3	Hook & Loop straps Width top/foot: 38-25cm Length: 56cm	1-Chamber Weight: 0,7kg	
	Arm Splint, angle VS56X01 VS56CX01 VS56CHX01 VS57CPX01 VS57CPHX01 Restraint straps: Width top/hand:	Red Stripes with Webbing CLICK with Webbing CLICK with Wipe-Clean Straps CLICK with Webbing CLICK with Wipe-Clean Straps	1-Chamber 1-Chamber 1-Chamber DOT-Chamber DOT-Chamber	
	Leg Splint VS58X02 VS58CX02 VS58CHX02 VS59X01 VS59CX01 VS59CHX01 VS59CPX01 VS59CPHX01 Restraint straps:	Red Stripes with Webbing CLICK with Webbing CLICK with Wipe-Clean Straps Red Stripes with Webbing CLICK with Webbing CLICK with Wipe-Clean Straps CLICK with Wipe-Clean Straps CLICK with Wipe-Clean Straps CLICK with Wipe-Clean Straps 4 Width top/foot: 38-25cm Length: 56	1-Chamber 1-Chamber 1-Chamber Multi-Chamber-Sack Multi-Chamber-Sack Multi-Chamber DOT-Chamber DOT-Chamber Cm Weight: 0,7kg	

RedVac Pediatric Systems

for use on the ambulance gurney


Pediatric Vacuum-Mattresses



VKF04X01 Cleats with Webbing Multi-Chamber-Sack
VKF04CX01 CLICK with Webbing Multi-Chamber-Sack
VKF04CHX01 CLICK with Wipe-Clean Straps Multi-Chamber-Sack
VKF04CPX01 CLICK with Webbing DOT-Chamber
VKF04CPHX01 CLICK with Wipe-Clean Straps DOT-Chamber

Restraint straps: 5 Securing straps: 5

Width Head/Middle/Foot: 53-61-53 cm Length: 140cm Weight: 3,5 kg



VKF06X02 Cleats with Webbing
VKF06CX02 CLICK with Webbing
VKF06CHX02 CLICK with Wipe-Clean Straps
VKF06CPX02 CLICK with Webbing
VKF06CPHX02 CLICK with Wipe-Clean Straps

CLICK with Wipe-Clean Straps

CLICK with Webbing

CLICK with Wipe-Clean Straps

Multi-Chamber-Sack

DOT-Chamber

DOT-Chamber

Multi-Chamber-Sack

Multi-Chamber-Sack

Carry handles: 4 Restraint straps: 3

Width Top/Hand: 85-65 cm Length: 115cm Weight: 2,9 kg

RedVac JUNIOR

Pediatric Restraint System for use with or without vacuum mattress





JUNIOR is a restraint system for safe transportation of children in ambulances. It is compatible with most ambulance gurneys. The system is suitable for new-borns and children up to around 7 years with body length 130 cm and weight of 35 kg.

Junior	KF01X07	
	Width head: 39 cm Width foot: 20 cm Length: 74 cm Weight: 1,3 kg Restraint straps: 4 Securing straps: 6	
Vacuum mattresses for JUNIOR	VKF01PX02 Hook & Loop straps DOT-Chamber VKF01CPX02 CLICK with Webbing DOT-Chamber VKF01CPHX02 CLICK with Wipe-Clean straps DOT-Chamber Width: 45 cm Length: 80 cm Weight: 0,75 kg Restraint straps: 3	
RAAA	VKF02PX01 Hook & Loop straps DOT-Chamber VKF02CPX01 CLICK with Webbing DOT-Chamber VKF02CPHX01 CLICK with Wipe-Clean straps DOT-Chamber Width Head/Middle/Foot: 45/62/62 cm Length: 138 cm Weight: 1,2 kg Restraint straps: 5	

RedVac*

Accessories for Vacuum Products

Can only be ordered together with the product

Antiskid systems

AS01

Three (3) stars with six cleats each



AS02 Three (3) rubber bands



Integrated Head Immobiliser

HI01







Integrated Leg Separator

LS01





Customised labels

Property indicator label



Label sizes:

8cm x 5cm

17cm x 7cm

Α4

Repair Kit

9022 Repair Kit with self-adhesive patches





V/		
Vacuum Pump	s & Carry	y Bags
	HP01 One-Hand pump,	plastic
		Width: 14cm Weight: 0,6kg
	HP03 Hand pump, alum Length: 30cm	ninium Width: 4cm Weight: 0,4kg
	HP04 Hand pump, plas	•
	FP01 Foot pump, plast Length: 43cm	ic Width: 11cm Weight: 0,75kg
a p	HFP02 Hand- and Foot p Length: 45cm	oump, plastic. Width: 6cm Weight: 0,95kg
	FP04 Foot pump, alum Length: 54cm	inium Width: 6cm Weight: 1,4kg
Carry bag for RedVac Junior	TS04	Nylon fabric
A	TS04C B/R	Colour: blue Cordura fabric
	TS04V B/R/G/Y	Colour: blue, red Wipe-off foil Colour: blue, red, green, yellow
	Length-Width-He	ight: 40-20-10cm Weight: 0,25kg
Carry bag for splint set	TS06	Nylon fabric Colour: blue
	TS06C B/R	Cordura fabric Colour: blue, red
	TS06V B/R/G/Y	Wipe-off foil Colour: blue, red, green, yellow
	Length-Width-He	ight: 65-45-12cm Weight: 0,4kg
Carry bag for mattress with backpack straps	TM01	Nylon fabric Colour: blue
	TM01C B/R	Cordura fabric Colour: blue, red
	TM01V B/R/G/Y	Wipe-off foil Colour: blue, red, green, yellow
	for mattress with Length-Width-He	width less than 90cm ight: 80-45-25cm Weight: 0,7kg
	TM02	Nylon fabric Colour: blue
	TM02C B/R	Cordura fabric Colour: blue, red
	TM02V B/R/G/Y	Wipe-off foil Colour: blue, red, green, yellow
		width of more than 90cm ight: 100-45-30cm Weight: 1kg

TRS Thermo Rescue Sack



for use with a vacuum mattress and/or spine boards or rescue stretchers

- Water repellent and rugged material provides good protection without additional blanket
- Designed for easy access to the patient
- 10 carry handles
- Under surface with antiskid stripes
- Zipper and Hook-and-Loop closing for fast handling and easy access to the patient
- Optional property indicator label














Auch sehr große Patienten können aufgenommen werden


Specification		
	TRSX01	
	Material: Nylon	
	Carry Handles: 10 Width Head/Middle/Foot: 45cm/70cm/40cm Length: 200cm Weight: 2,5kg	
	TRSX02	
	Material: Nylon	
	Carry Handles: 10	
	Width Head/Middle/Foot: 45cm/70cm/40cm Length: 220cm Weight: 2,9kg	
	TRSX03	
	Material: CORDURA quilted	
	Carry Handles: 10 Width Head/Middle/Foot: 45cm/70cm/40cm Length: 220cm Weight: 2,9kg	

BRS Helicopter Rescue Sack



for use with a vacuum mattress and/or spine boards or rescue stretchers





- Antiskid under surface to prevent slipping on snow
- Hook and loop fastener and zipper for easy access to the patient and to avoid wind-sack phenomenon and rotation
- Stable head protection
- Lifting strap system is adjustable to achieve horizontal position of the patient
- Patient restraint straps with metal buckles
- Pocket for vacuum pump
- Plated-through holes for free drainage of water

Anti-rotation rings

CE labelled and certified for helicopter lifting




BRS03

All straps are sewn onto the rescue bag and are stowed in easy to open pockets



The detachable lifting bridle is attached to the carry handles with karabiner



Integrated bag for the vacuum pump





Undersurface with antiskid stripes



Specification



BRS03 BRS04

Carry Handles: 10

Patient restraint straps: 5

Width, Head-Middle-Foot: 40-70-45cm

Length: 220 cm Weight: 4 kg

TBRS03

Carry bag

Width-Length-Height: 44-80-33 cm

Weight: 1 kg



Quality management in accordance with EN ISO 9001 - EN ISO 13485





KOHLBRAT & BUNZ GMBH

A-5550 Radstadt Austria

T: +43 6452 7193 0

F: +43 6452 7193 51

office@kohlbrat-bunz.com www.kohlbrat-bunz.com

www.redvac.com